Short Manual / Kurzanleitung / Bref manuel d'instructions

peqSTAR

Thermal cyclers / Thermocycler / Thermocycleur

peqSTAR XS 32 well thermal cycler

peqSTAR 2X double block thermocycler

peqSTAR 96X, 384X and in situ X thermocycler /

peqSTAR 32 Well-Thermocycler

peqSTAR 2X Doppelblock-Thermocycler

peqSTAR 96X, 384X und in situ X Thermocycler

CONTENTS

1	Safety instructions	2
2	License Notice	2
3	System Overview	3
4	Installation	3
4.1	Content of delivery	3
4.2	Installing the thermal cycler	3
4.3	Operation of the thermal cycler	3
5	Cycler Operation	4
5.1	Start screen = 'Run'	4
5.2	Running stored programs ('Run'):	4
5.3	Incubation mode ('Run')	5
5.4	Create or edit PCR programs ('Programs')	5
5.5	Further dialogues: Diagnostics, GLPs, System	9
5.6	Further questions...?	10
6	Technical Service and Ordering Information	11
7	Technical data	12
8	Maintenance and repair of the Thermal cycler	14
8.1	Cleaning, maintenance and repair of the thermal cycler	14
8.2	Over temp test	14
8.3	Temperature check-up	14
8.4	Calibration	14
9	Warranty	15
10	Disposal Instruction	15

A more detailed instruction manual in English or German language is available for download on our website www.peqlab.com.

1 SAFETY INSTRUCTIONS

Before the first use of the thermal cycler please read this instruction manual entirely.

Special notice should be taken of the following

Caution of dangerous voltage

Please ensure the voltage indicated on the device exactly matches your local electrical supply

Caution of dangerous explosive material

Explosive or reactive material mustn't be heated in the thermal cycler.

Caution - liquids

Ensure that no liquids can enter the device.

Caution - hot surfaces

The thermal plate and the heated lid quickly reach temperatures above 50 °C. There is danger of burning! Keep the lid closed until the temperatures reaches 30 °C or less.

Only use materials (tubes, tube caps, plates and sealing films) which are heat resistant to 120 °C

Caution – environment

The ventilation of the device must not be covered.

2 LICENSE NOTICE

This instrument is licensed for research and development and for uses other than human in vitro diagnostics under one or more of the following patents of Applied Biosystems Corporation: U.S. Patent Nos. 5,656,493, 5,038,852, 5,333,675, 5,475,610 (claims 1-159 and 164-166), 6,703,236 (claims 7-10) and 7,238,517 or corresponding claims in their non-U.S. counterparts. No right is conveyed expressly, by implication or by estoppel under any other patent of Applied Biosystems, including but not limited to U.S. Patent No. 6,814,934 and its non-U.S. counterparts, which describe and claim thermal cyclers capable of real-time detection.

3 SYSTEM OVERVIEW

The peqSTAR thermal cyclers are the ideal PCR devices for all molecular biological and biochemical laboratories working in the field of basic research or routine diagnostic. Eight or sixteen 250 °C HTR (High Temperature Range) Peltier elements allow quick changes in temperature up to 5 °C/sec (max) (3°C/sec at peqSTAR XS). The individual monitoring and controlling of each of the eight Peltier elements equipped with 'Long Life Technology' ensure the tightest correlation of target and actual temperature across the entire thermal plate leading to highest reproducibility of the results. Each block is supplied with a lid heating to avoid condensation (except for in Situ blocks).

The peqSTAR thermal cycler is also equipped with a touch-sensitive, graphic colored TFT display of the newest generation, allowing operation of the device by simply touching symbols and graphic elements on the screen. This leads to easy and intuitive handling and programming of the device. Additionally a mouse can be connected via USB. Further-more an Ethernet port is available to connect the device with the Microsoft Windows® or Linux network. It is therefore possible to use an external server for archiving programs and GLP reports which will be available globally in the network (e.g. for other thermal cyclers).

It is also possible to control the thermal cycler by a PC integrated to the network or wireless ZigBee connection. The provided software offers the same functions as the internal thermal cycler software, except for some system settings. Transfer of data can be done by using commercially available USB memory sticks enabling a nearly unlimited saving space for programs and GLP reports, which can be printed using a printer connected to the thermal cycler via the network.

Notice: The USB ports support only standard cable mice and standard USB sticks.
The maximum acceptable length of the connecting cable of the port is 3 m!

4 INSTALLATION

4.1 Content of delivery

- 1 Thermal cycler
- 1 Power cord
- 1 Short manual

(A detailed instruction manual in English or German is available for download on our website www.peqlab.com.)

- 1 CD PC-software 'peqSTAR Manager'

Please check delivery for completeness and transport damage upon arrival. If any transport damage is noticed, contact your distributor or manufacturer immediately.

NOTICE: IF ANY DAMAGE OF THE DEVICE SHOULD BE FOUND DO NOT USE THE DEVICE!

4.2 Installing the thermal cycler

Check the packaging for any transport damage. Remove the packaging and place the thermal cycler onto a solid surface. The device should not be exposed to direct sunlight. There should be enough space to make sure that the ventilator grill on the front and back are not covered and sufficient ventilation is provided. There should be a distance of at least 25 cm to the wall or the next device. Two thermal cyclers should not be placed back-to-back. Room temperature between 4 °C und 25 °C is advised.

All original packaging should be stored in case a need arises to return the instrument to the supplier.

4.3 Operation of the thermal cycler

Before starting up the thermal cycler for the first time make sure that the voltage range of the device is suitable to the electricity supply. Plug in the mains plug of the thermal cycler to the electricity supply. The outlet and the power switch are located at the back of the device. The serial number is written onto a typed label at the bottom of the device.

The device should only be used by authorized staff. To ensure this, a user administration with appointing rights for single users is available in the device.


5 CYCLER OPERATION

The operation of the thermal cycler can be run by the internal touch screen, either using the touch sensitive surface of the display or by using a mouse connected to a free USB port.

If the user has to enter numbers or text the corresponding field has to be touched on the touch screen or clicked with the mouse. An appropriate keyboard will appear on the screen according to the expected entry.

5.1 Start screen = 'Run'



The last ran PCR program is displayed in the 'Run' dialogue and can directly be restarted with the button 'Start PCR'.

Note: At peqSTAR 2X double block thermal cyclers the desired block can be chosen in the dialogues 'L' (left) and 'R' (right). For peqSTAR XS, 96X, 384X and in Situ X systems this option is not available.


5.2 Running stored programs ('Run'):

1. Open the selection dialogue in the 'Run' dialogue with the button .
2. Select the desired program in the selection dialogue and load to the 'Run' dialogue with the button 'Open'.
3. The loaded program is started by using the  button.



5.3 Incubation mode ('Run')

Use the button  in the 'Run' dialogue to set the device into the incubation mode. The temperature for the thermal plate and lid can be set manually. The clock shows the time passed by in the incubation mode.

5.4 Create or edit PCR programs ('Programs')

1. To create or edit PCR programs go to the 'Programs' dialogue.
2. By touching the 'Programs' dialogue the first time during the session the selection dialogue shows up where a new program can be set over 'New' or an existing program can be opened over 'Open' to edit.
3. On the left side of the program editor all available commands are displayed, on the right side of the editor the list of commands (Steps) the program is composed of is shown. After selecting a program command out of the list in the left window it can be inserted into the desired place via the  button.
4. Via double click or the 'Edit' button  parameters of single steps can be changed later.
5. Selected program steps can be deleted from the program by using the 'Delete' button .
6. Use the 'Save' button  to save the program.
7. Save and start programs in one step by simply using the 'Save&Run' button .
8. For even faster programming the peqSTAR thermocycler comes with a selection of template programs, which could easily be edited and saved under a new program name.

Adding and removing program  Program name  Program window steps

Available commands:


	Lock Lid Locks the lid avoid accidental opening.		Cycle Combines Open Loop, Close Loop, Temperature, Ramp and Gradient.
	Unlock Lid Unlocks the lid in case that the lid has been locked in a previous program step.		Pause Interrupts the program temporarily at the definable temperature.
	Heat Lid Turns the lid heating on/off.		Store Holds the definable temperature, turns off the heated lid automatically.
	Temperature Temperature step for definable time, Temperature and Time Incr./Decr. function.		Start Cycle Manual loop beginning.
	Gradient Realizes different temperatures at the same time for protocol optimization (see 'FlexGradient' below for more information).		Close Cycle Manual loop ending.


Sample program

Create a new program with 'New'. Enter a program name and confirm with **OK** to access the programming functions.


01	 Heat Lid	1 Heat Lid to 110.0°C	Heated Lid on, temperature 110 °C A 'Heat Lid' step is inserted in each program by default as first step. This step can be edited or deleted later. The standard settings could be changed under 'System' => 'Firmware' => 'Settings'. For 'Store' or 'Temp' steps < 30 °C the lid heating is automatically switched off.
----	---	-----------------------	---

To enter another program step select the desired command out of the list of the available command and add it to the program window with the '+' arrow. Thereupon the corresponding parameters will be asked.

02	 Temp	2 Temp. 95.0°C for 2' 0"	Denaturation at 95 °C for 2 minutes
03		3 Start Cycle, 30x	Open the 'Cycle' function
-		4 Temp. 95.0°C for 30"	→ Enter number of cycles – 30 cycles
07	Cycle	5 Temp. 65.0°C (1.0°C/s) for 30" (-... 6 Temp. 72.0°C for 1' 0" (+0.0°C,+... 7 Close Cycle	→ 3 x 'Temp' step → Terminate 'Cycle' function with 'Close Cycle'
04	 Temp	4 Temp. 95.0°C for 30"	Denaturation at 95 °C for 30 seconds

05	 Temp	<p>5 Temp. 65.0°C (1.0°C/s) for 30" (...)</p>	<p>Annealing at 65°C in the 1st cycle for 30 seconds 'Touchdown'-PCR via Temperature-'Decrement'. The annealing temperature of 65 °C in the 1st cycle is reduced by 0.4 °C each cycle. This leads to high specificity in the first cycles and high amplification efficiency in the later cycles.</p> <p>The reduction of the cooling rate at this step of the program is not necessary, but should only be an example. It might make sense e.g. for quantitative annealing of complementary oligo nucleotide or at annealing steps in the first cycles, when primer annealing should be very efficient because of low template amount (e.g. after cDNA synthesis)</p>
06	 Temp	<p>6 Temp. 72.0°C for 1' 0" (+0.0°C,+...</p>	<p>Elongation at 72 °C for 1 minute in the 1st cycle Via Time -'Increment' the elongation time is prolonged by 5 s each cycle. The incremental prolongation of the elongation time is useful with high numbers of cycles (>30), when the Polymerase activity declines.</p>
07	 Close Cycle	<p>7 Close Loop</p>	<p>Terminate 'Cycle' function with 'Close Cycle'.</p>

Alternative programming of steps 03 – 07:


03	 Start Cycle	<p>3 Start Cycle, 30x</p>	<p>Open program cycle – 30 cycles.</p>
04	 Temp	<p>4 Temp. 95.0°C for 30"</p>	<p>Denaturation at 95 °C for 30 seconds</p>
05	 Temp	<p>5 Temp. 65.0°C (1.0°C/s) for 30" (...)</p>	<p>Annealing at 65°C in the 1st cycle for 30 seconds 'Touchdown'-PCR via Temperature-'Decrement'. The annealing temperature of 65 °C in the 1st cycle is reduced by 0.4 °C each cycle. This leads to high specificity in the first cycles and high amplification efficiency in the later cycles.</p> <p>The reduction of the cooling rate at this step of the program is not necessary, but should only be an example. It might make sense e.g. for quantitative annealing of complementary oligo nucleotide or at annealing steps in the first cycles, when primer annealing should be very efficient because of low template amount (e.g. after cDNA synthesis)</p>
06	 Temp	<p>6 Temp. 72.0°C for 1' 0" (+0.0°C,+...</p>	<p>Elongation at 72 °C for 1 minute in the 1st cycle Via Time -'Increment' the elongation time is prolonged by 5 s each cycle. The incremental prolongation of the elongation time is useful with high numbers of cycles (>30), when the Polymerase activity declines.</p>
07	 Close Cycle	<p>7 Close Cycle</p>	<p>Terminate program loop with 'Close Cycle'.</p>
08	 Temp	<p>8 Temp. 72.0°C for 5' 0"</p>	<p>Final elongation at 72 °C for 5 minutes The final elongation assures the complete synthesis of all amplificats.</p>
09	 Store	<p>9 Store forever at 8.0°C</p>	<p>8 °C forever, lid heat is turned off By the 'Store' step the lid heat is turned off automatically. This avoids that lid heat and Peltier elements of the block work against each other. Especially at longer incubation times at low temperature this increases the life time of the Peltier elements and reduces the power consumption.</p>

FlexGradient Technology™


The temperature of each of the eight rows of the block is individually controlled! For the first time you can choose either the perfectly linear temperature gradient – ideally suited for PCR optimizations, or independent lane control – ideally suited for the use of different primer pairs in one PCR run.

First, the temperature range is defined by the buttons 'Center' and 'Gradient ±'. The standard gradient is always realized from the top to the bottom over the 8 rows of the block. Moreover with the peqSTAR 96X you have the opportunity to set a 16-gradient. In this case different temperatures are realized in the left and right part of the thermoblock.

As in this case you will also have a gradient in the lanes from left to right, the temperatures displayed in the software are only realized in the rows 3 and 10.

By selecting one of the temperature buttons **57.9** different options will be available: When the checkbox linear was deactivated before, the single temperatures of each row can be adjusted within a special range using the button 'Adjust Temperature'. Moreover the selected temperature can be used as new center value or the gradient step can be changed in a temperature step using the selected temperature.

Program options:

When setting up a new program, when saving or during programming (via the Info button) the program options can be set. The following options are available for the user.

1. Power Fail Denaturation:

If a power failure had happened while a program was running the thermal cycler restarts the program at the point it was stopped as soon as the power returns.

If the option 'Power Fail Denaturation' is activated the thermal block will heat to the stated temperature for the set time before continuing with the program.

2. Global Program Ramp:

Using this function a global ramp rate in °C/s can be programmed. This global ramp rate will be used in all temperature and gradient steps of the PCR program, but could be changed in a single temperature or gradient step manually by entering the desired ramp value at the step.

3. Gradient Control

Using this option you can decide whether in a gradient step the plateaus (Plateaus synchronized) or the Ramps (Ramps synchronized) should be synchronized for all rows.

4. Tube Control

Via tube control you have the possibility to adjust temperature regulation for a temperature step according to sample volume. If necessary, the temperature regulation will perform short term deliberate overshoots and undershoots depending on the entered sample volume during PCR run.


Especially if combined with the emulation mode this function can be helpful if the emulated cycler used the tube control function before, too.

5. Emulation


Using the emulation function the thermal cycler will perform PCR like the chosen cycler, emulating its performance. This function could be useful if a PCR program should be transferred to the peqSTAR 2X from another type and model of thermal cycler.

5.5 Further dialogues: Diagnostics, GLPs, System


Diagnostics: Monitoring the actual temperatures of each of the six Peltier elements during a run. In suspicion case of an instrument malfunction have a look at the Peltiers' temperature profile.

GLPs: Displaying and organizing GLP (Good Laboratory Practice) reports created by the thermal cycler. This option is mainly for labs who need to save GLP reports for documentation purpose.

System: Configuration of the system and changing of the general settings of the device.

Elements of the system dialogue:

	Time Setting up date and time.		Hide/Show Cursor Hide or show the cursor on the screen.
	User Call Setting of acoustic signals or email notification (Pause, Store, Program End)		User Management Define user and user groups (and their access rights). The default password for the user 'Administrator' is: ' admin '. Via 'Logout' button the current user can be logged off. Hereby the user 'Guest' is automatically logged in. This user only has minimal user rights. These settings can be changed by the administrator if desired.
	Network settings Network settings for network integration of the peqSTAR.		
	Hardware Display hardware versions of Control Unit, Power Unit and Thermo block.		Firmware Display, administrate and update software (Control/Power Unit, Thermo block)
	Printer Set up a printer path for an existing network printer.		Display Settings Change settings of the screen display.
	Protocol Internal communication of the thermal cyclers for failure diagnosis		Service Service function (only for authorised technical Service)
	LAN Control Functional button for PC and Master/Slave control		USB Recovery Manual mounting of USB sticks that are not automatically detected.

5.6 Further questions...?

Direct help :

Use this button to activate (deactivate) the direct help mode. If the mouse pointer is visible a question mark appears alongside. By touching a function button or other dialogue element the corresponding help window will appear.

A more detailed instruction manual in English or German is available for download on our website www.peqlab.com.

6 TECHNICAL SERVICE AND ORDERING INFORMATION

Should you have any questions about the peqSTAR thermocycler, please visit www.peqlab.com to find the respective contact person.

This short manual is delivered with the following products:

Product		Cat. No.
peqSTAR XS	with 32 well Universal block with standard lid for 32 x 0.2 ml tubes or 16 x 0.5 ml tubes with flat caps.	95-03002
peqSTAR 2X Universal	with 2 x 48 well Universal blocks with standard lids for 48 x 0.2 ml tubes, 48 well PCR plates or 24 x 0.5 ml tubes with flat caps each	95-07002
peqSTAR 2X Universal Gradient	with 2 x 48 well Universal Gradient blocks and standard lids for 48 x 0.2 ml tubes, 48 well PCR plates or 24 x 0.5 ml tubes with flat caps each	95-08002
peqSTAR 96X Universal	with 96 well Universal block and standard lid for 96 x 0.2 ml tubes, 96 well PCR plates or 48 x 0.5 ml tubes with flat caps	95-05002
peqSTAR 96X Universal Gradient	with 96 well Universal Gradient block and standard lid for 96 x 0.2 ml tubes, 96 well PCR plates or 48 x 0.5 ml tubes with flat caps	95-06002
peqSTAR 96X HPL	with 96 well block and High Pressure Lid for 96 x 0.2 ml tubes or 96 well PCR plates	95-05002-HPL
peqSTAR 96X HPL Gradient	with 96 well Gradient block and High Pressure Lid for 96 x 0.2 ml tubes or 96 well PCR plates	95-06002-HPL
peqSTAR 384X HPL	with 384 well block and High Pressure Lid for 384 well PCR plates	95-050384-HPL
peqSTAR <i>in situ</i> X	with <i>in situ</i> block for up to four microscope slides and integrated buffer reservoir	95-050FLAT

7 TECHNICAL DATA

General characteristics

- 8, 2 x 8 or 16 Peltier elements with Long-Life-technology and 8, 2 x 8 or 16 control circuits with Pt 1000 temperature sensors, respectively (at peqSTAR XS, 2X and 96X)
- Maximal heating and cooling rate: 5 °C/s (3°C/s for peqSTAR XS)
- Block uniformity (at 72 °C): ± 0.20 °C
- Temperature range thermoblock: 4 to 105 °C
- Control accuracy thermoblock: ± 0.1 °C
- Variable ramping: 0.1 to 3.0 °C/s
- Increment/Decrement time: 0:01 to 9:59 minutes
- Increment/Decrement temperature: 0.1 to 9.9 °C
- Lid heating temperature range: 40 to 120 °C
- Programmable Lid locking mechanism (except for peqSTAR XS)
- 4 x USB, 1 x Ethernet interfaces (MS Windows® or Linux), maximal length of power cord 3 m!
Note: USB interfaces only support standard USB sticks, cord mouse and keyboards!
- Dimensions (B x H x D): 30 x 28 x 38 cm
- Dimensions peqSTAR XS (B x H x D): 27 x 19 x 33 cm
- Weight incl. block: 13.3 kg (5.2 kg for peqSTAR XS)
- Electrical power supply: 90 – 240 V AC, 50 – 60 Hz, 850 VA (150 VA for peqSTAR XS)
- Fuse: 10 A time lag
- Degree of pollution: 2
- Ambient temperature: 10 °C to 30 °C
- Maximum relative humidity: 70 %
- Maximum altitude for the use of this system: up to 2000 m above sea level
- For indoor use only!
- Noise level: < 37 dB (A)

Gradient feature (optional)

- Maximal gradient range over 8 or 16 columns, respectively: 30 °C (± 15 °C)
- Temperature range Gradient: 35 to 105 °C
- Gradient accuracy: ± 0.1 °C
- Gradient function can be upgraded even after purchase and on-site
- FlexGradient™ technology: ideal, linear Gradient or independent lane control
- 'Gradient Control': synchronized plateau times or ramp rates

User interface/functions

- Touch sensitive TFT-Display (VGA, Graphic, 65535 colours), operation via USB mouse possible
- Direct help function, Tm- and elongation time calculator
- Internal flash memory for 500000 typical PCR programs in free configurable folders/subfolders
- Unlimited number of programs via network PC or USB memory stick
- Free computer Software for the generation of PCR programs on the computer
- Unlimited number of steps/program
- Password protected user accounts with variable access rights
- GLP reports for complete recording of all runs
- Quickstart function of the last used program
- Automatic restart after power failure, Power-Fail-Denaturation, Instant Incubation
- Remote control of the system by PC software
- Free choice of mp3 signal sounds
- User notification by email
- Master/Slave control (wireless option)

Block options (partly available as interchangeable blocks)

peqSTAR XS:

32 well Universal block with: lid heating and automatic height adaption for 32 x 0.2 ml tubes or 16 x 0.5 ml tubes with flat caps.

peqSTAR 2X/2X Gradient:

Two 48 well Universal blocks with: lid heating and automatic height adaption for 48 x 0.2 ml tubes, 48 well PCR plates or 24 x 0.5 ml tubes with flat caps.

peqSTAR 96X/96X Gradient:

96 well Universal block with: lid heating and automatic height adaption for 96 x 0.2 ml tubes, 96 well PCR plates or 48 x 0.5 ml tubes with flat caps.

peqSTAR 96X HPL/96X HPL Gradient:

96 well block with High Pressure lid (HPL, 100 – 250 N) for the safe sealing of 96 well PCR plates

peqSTAR 384X HPL:

384 well block with High Pressure lid (HPL, 100 – 250 N) for the safe sealing of 384 Well PCR plates

peqSTAR In situ X:

In situ block with integrated buffer reservoir, for up to 4 microscope slides

8 MAINTENANCE AND REPAIR OF THE THERMAL CYCLER

8.1 Cleaning, maintenance and repair of the thermal cyclers

For safety reasons the device must be switched off and the mains plug pulled before cleaning is carried out!

The surface of the thermal cyclers can be wiped off with a damp cloth or ethanol. Do not use aggressive or scouring cleaners or organic solvents for cleaning. The device should be protected from aggressive chemicals. Make sure that no liquid penetrates into the interior of the device.

If sample liquids penetrate into the drill holes of the thermal cyclers clean them immediately with mild soap followed by distilled water. You can use mild disinfectants for decontamination.

Maintenance (Overtemp test and temperature check-up, see 8.2 and 8.3) should be carried out regularly to detect possible technical failures.

The electric fuses are placed at the back of the device between the power switch and the mains plug.

After pulling off the main plugs and pushing upwards a safety catch the fuses can be changed. Use only fuses with correct values (indications about the fuse type are placed at the back of the device).

The device may only be opened by a qualified specialist. Unauthorized work on the device voids the warranty. All kinds of repairs may only be carried out by authorized persons using original replacement parts. The replacement of single components of the device should not be done by the user but exclusively by authorized specialists.

Therefore the device must be sent to the manufacturer. The surface of the device should be decontaminated with a lint free cloth soaked with 70 % ethanol before sending. The drill holes of the thermoblock shall be cleaned with a 1.5 % hypochlorite solution with the help of a cotton bud. Please fill in and sign the decontamination certificate and send it to the manufacturer together with the device.

8.2 Over temp test

A periodical test of the overtemperature protection (→ System → Hardware → Thermo Block → Overtemperature Test) is recommended.

8.3 Temperature check-up

Maintenance of the temperature at regular intervals (e.g. annually) is recommended. The thermometer for use should be equipped with wired sensing devices (e.g. Pt1000) and should have an accuracy of measurements of 1/10 class B ($\Delta T = \pm 0.03 \text{ }^\circ\text{C}$) in the range of 0 to 100 °C. While processing the temperature control the room temperature should be 25 °C at maximum.

8.4 Calibration

If there is a reasonable suspicion that the device has deviated from the calibration range, because the variations of the temperature check-up are considerably out of the acceptable tolerance, a new calibration must be carried out by the manufacturer. Please contact your local distributor.

9 WARRANTY

The period of warranty for the product at hand is 24 month from date of purchase. The warranty does not apply for defects caused by incorrect, improper or inappropriate use. The customer has to check up on obvious defects promptly upon receipt of the device and send a written report. Otherwise the claim for warranty expires. The manufacturer reserves the right to carry out changes to the thermal cycler or to the technical documentation of the thermal cycler on the basis of product development or improvement without prior notice.

10 DISPOSAL INSTRUCTION

This product is subject to the WEEE disposal regulations and cannot be disposed of with regular waste.

Please contact VWR International GmbH or your local distributor for disposal of peqSTAR thermal cyclers. Waste electronic equipment is not useless waste. Environmentally friendly disposal can help to retrieve valuable resources. Thereby you contribute to the protection of the environment and human health.

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This instrument is licensed for research and development and for uses other than human in vitro diagnostics under one or more of the following patents of Applera Corporation: U.S. Patent Nos. 5,656,493, 5,038,852, 5,333,675, 5,475,610 (claims 1-159 and 164-166), 6,703,236 (claims 7-10) and 7,238,517 or corresponding claims in their non-U.S. counterparts. No right is conveyed expressly, by implication or by estoppel under any other patent of Applera, including but not limited to U.S. Patent No. 6,814,934 and its non-U.S. counterparts, which describe and claim thermal cyclers capable of real-time detection.

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