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Thermo Shaker TS-100 / Thermo-Inkubationsmischer TS-100

CONTENTS

1.	WARRANTY	1
2.	SAFETY PRECAUTIONS	1
3.	GENERAL INFORMATION	2
4.	INSTALLATION	2
5.	OPERATION	3
6.	TECHNICAL DATA	4
1.1	Temperature specifications	4
1.2	General specifications	4
7.	CLEANING AND MAINTANCE	5
8.	TECHNICAL SUPPORT	5

1. WARRANTY

PEQLAB guarantees that the system you have received has been thoroughly tested and meets its published specification.

However, immediately upon arrival, please check carefully that the shipment is complete and has not been damaged in transit. For missing parts or to report any kind of damage, please contact PEQLAB (see 'TECHNICAL SUPPORT'). Please retain all packaging materials until the delivery has been completely checked since this will speed up the return of goods if required and reduce environmental impact. Any form of returns, replacements or credit notes must be agreed in advance by PEQLAB.

For the Thermo Shaker TS-100 PEQLAB guarantees a warranty period of 24 months if the products have been used solely according to the instruction manual and if not agreed differently. No liability is accepted for loss or damage arising from incorrect use. PEQLAB's liability is limited to the repair or replacement of the unit or refund of the purchase price, at PEQLAB's discretion. PEQLAB is not liable for any consequential damages.

PEQLAB reserves the right to alter the technical specifications without prior notice. This will enable us to implement developments as soon as they arise.

2. SAFETY PRECAUTIONS

The following symbols mean:

Caution: Read these operating instructions fully before use and pay particular attention to sections containing this symbol!

Caution: Surfaces can become hot during use!

- Use only as specified by the operating instructions, or the intrinsic protection may be impaired.
- Connect only to a power supply with a voltage corresponding to that on the serial number label.
- Do not use other power supply units than recommended by the manufacturer.
- Ensure that the mains switch and isolating device (power supply connector) are easily accessible during use.
- Connect only to a power supply which provides a safety earth (ground) terminal.
- Do not operate the unit outside the laboratory premises.
- Do not leave the operating unit unattended.
- Do not operate the unit in premises with aggressive or explosive chemical mixtures.
- Before moving, disconnect at the power supply socket.
- Do not impede the platform motion during operation.
- Use only standard qualitative tubes.
- If liquid is spilt inside the unit, disconnect it from the power supply and have it checked by a competent person.
- Before using any cleaning or decontamination method except those recommended by the manufacturer, user should check with the manufacturer that the proposed method will not damage the equipment.
- The unit should be saved from shocks and falling.

3. GENERAL INFORMATION

TS-100 is a well-established popular model of the thermal shaker for microtest tubes. It is the ideal instrument for intensive mixing of samples in the regulated temperature conditions.

Temperature control and mixing can be used both simultaneously and independently i.e. the unit implements three products in one: 1) a shaker, 2) a thermostat and 3) a thermal shaker.

TS-100 is applicable in:

- DNA analysis for DNA extraction and further sample preparation
- Biochemistry for studies of enzymatic reactions and processes
- Metabolite extraction from cell materials, polysaccharides, lipids and other cell components

Thermo-Shakers provides:

- Soft, but intensive rotational stirring of the samples
- Smooth regulation, stabilization and indication of rotation speed
- Even amplitude throughout the shaker platform
- Setting up and indication of the required working time
- Automatic stop of the platform motion after the set time expires
- The set parameters preservation in memory
- Indication of the current working time
- Set up and indication of the set and actual temperature on the block

4. INSTALLATION

Unpacking: Remove packing materials carefully and retain for future shipment or storage of the unit.

Thermo-Shaker set includes:

Thermo-Shaker TS-100, 2 spare rubber belts, power cord, AC/DC Adapter, Operating Manual

Place the unit upon even horizontal non-flammable surface away from any flammable materials. To provide optimum ventilation ensure clearance around the device: 20 cm on each side. Ensure clearance on the underside, do not place any object under unit. Plug the external power supply unit (AC/DC adapter supplied with the shaker) into the 12 V socket at the rear side of the Thermo-Shaker.

Block installation (if block is not installed at the moment of delivery):

Connect the plug (fig. 1/2) to the contact terminal on the underside of the block. Align the block so that the warning label is facing the front of the unit and secure with the four knurled screws.


Figure:1

Changing blocks:

- Disconnect the external power supply unit from the mains.
- Remove the four knurled screws (fig.1/1), disconnect the plug (fig.1/2).
- Select the new block.
- Connect the plug. Align the block so that the warning label is facing the front of the unit and secure with the four knurled screws.

5. OPERATION

Connect power supply unit to the mains and switch ON the power switch located on the rear panel of the unit.

The display will turn on with the upper line (set point) showing time, RPM and temperature set earlier and the lower line (actual point) showing current readings of the same parameters (STOP - time, 000 - RPM, thermoblock temperature °C which automatically starts rising according to the temperature set in the upper line). The time of temperature stabilisation depends on the initial temperature, but does not exceed 15 - 20 min, if the set temperature is 37.0 °C.

How to set the necessary parameters:

Use the readings in the upper line of the display (set point), while setting the necessary parameters.


Figure: 2

Reaction time (TIME):

With the help of "▲" and "▼" keys (Fig. 2/1) set the required working time interval in hours and minutes (increment - 1 min). If the key is pressed for longer time the increment becomes bigger.

Shaking intensity (RPM):

With the help of "▲" and "▼" keys (Fig. 2/2) set the required shaking intensity in revolutions per minute (increment 10 RPM). If the key is pressed for longer time the increment becomes bigger.

Reaction temperature (T, °C):

With the help of "▲" and "▼" keys (Fig. 2/3) set the necessary temperature (increment 0.1 °C). If the button is pressed for longer time the increment becomes bigger.

For TS-100 it is possible to turn off heating of the thermoblock by setting the required temperature below 25 °C (the display will show OFF - T, °C - set point). In this mode TS-100 can be used in the cold rooms as a mixing device without thermoregulation. The set parameters can also be changed during operation.

Program execution:

After the thermal stabilisation of the Thermo-Shaker (when the set and current temperature readings become the same):

Insert tubes into the platform sockets. Press the RPM-RUN/STOP key (Fig. 2/4). The platform will start rotation and the timer indicator will start counting up the time interval (with 1 min precision).

Please note: If the rotation speed is set to zero, pressing RPM-RUN/STOP keys starts the timer but the platform does not move.

At the end of the program (after the set time elapses) the platform motion stops and the timer shows the flashing reading STOP accompanied by the repetitive sound signal until the RPM-RUN/STOP key is pressed.

If the working time is not set (or deleted) and the timer indicator in the upper line shows 00:00, pressing the RPM-RUN/STOP button cause the Thermo-Shaker to operate continuously until the RPM-RUN/STOP button is pressed again.

If required, there is possibility to restart the timer when it is running. Press the TIME-RUN/STOP key once (Fig. 2/5) to stop the timer. Press the TIME-RUN/STOP key again to restart the timer.

The platform motion can be stopped at any time by pressing the RPM RUN/STOP key. In this case the program realisation and the platform motion stops and timer is set back to zero switching into the STOP mode. Press the RPM- RUN/STOP key to repeat the operation with the same time and speed.

At the end of operation turn OFF the unit with power switch at the rear panel and disconnect the external power supply unit from the mains.

6. TECHNICAL DATA

6.1 TEMPERATURE SPECIFICATIONS

- | | |
|---|-----------------|
| ▪ Temperature setting range (°C) | + 25 to +100 °C |
| ▪ Temperature uniformity over the platform (+ 25 to +40 °C) | ± 0.1 °C |
| ▪ Temperature uniformity over the platform (+ 40 to +70 °C) | ± 0.2 °C |
| ▪ Temperature uniformity over the platform (+ 80 °C) | ± 0.4 °C |
| ▪ Temperature uniformity over the platform (+ 100 °C) | ± 0.5 °C |

6.2 General specifications

Thermo-Shaker is designed for operation in closed laboratory rooms at ambient temperature from + 5 °C to + 40 °C and maximum relative humidity 80 % for temperatures up to 31 °C decreasing linearly to 50 % relative humidity at 40 °C.

- | | |
|--|------------------------------------|
| ▪ Speed range | 250 -1400 RPM (increment 10 RPM) |
| ▪ Orbit | 2 mm |
| ▪ Independent timer with sound signal | 1 min - 96 hrs (increment 1 min) |
| ▪ Continuous operation time | max. 96 hours |
| ▪ Time of thermoblock heating (RT to +37 °C) | 15 min |
| ▪ Display | 16 x 2 signs, LCD |
| ▪ Power | External power supply 12 V, 4.16 A |
| ▪ Dimensions | 205 x 130 x 230 mm |

- Weight (with power supply) 4.0 kg
- Standard thermoblock capacity:
 - TS-SC-20 20 x 1.5 ml tubes
 - TS-SC-20N 20 x 2 ml tubes
 - TS-SC-18 20 x 0.5 ml + 12 x 1.5 ml tubes
 - TS-SC-16 20 x 0.2 ml + 12 x 1.5 ml tubes

7. CLEANING AND MAINTANCE

If service is required, switch off the unit and contact peqlab for repairs. It is allowed to perform the maintenance and all types of repair operations only to the authorised specialists or service-engineers who have completed special training.

Cleaning liquids that do not contain concentrate organic solvents, alkali or acid can be used for device cleaning. Standard ethanol (75 %) can be used for disinfection.

Rubber belt replacement (fig.3):

Remove 4 fixation screws on the shaker bottom and remove the bottom plate. Replace the rubber belt. Assemble the unit.


Figure 3

8. TECHNICAL SUPPORT

For technical questions and more detailed information on PEQLAB's products please visit www.peqlab.com to find the respective contact person.

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D PEQLAB Biotechnologie GmbH, 91052 Erlangen, Freecall (D): 0800 100 20 16, info@peqlab.de, www.peqlab.de
AT PEQLAB Biotechnologie GmbH, 6404 Polling, Tel: +43 (0) 5238 84 169, info@peqlab.at, www.peqlab.at
UK PEQLAB Ltd., Southampton SO31 7ZN, Freephone (UK): 0808 202 1302, info@peqlab.co.uk, www.peqlab.co.uk
USA PEQLAB LLC, Wilmington, DE 19810, Toll-Free (US): 877 737 5220, info@peqlab.us, www.peqlab.us