IKA-ULTRA-TURRAX® T 25 basic

<i>BETRIEBSANLEITUNG</i>	<i>D</i>	<i>5</i>
<i>OPERATING INSTRUCTIONS</i>	<i>GB</i>	<i>17</i>
<i>MODE D'EMPLOI</i>	<i>F</i>	<i>29</i>
<i>INSTRUCCIONES DE MANEJO</i>	<i>E</i>	<i>41</i>
<i>HANDLEIDING</i>	<i>NL</i>	<i>45</i>
<i>ISTRUZIONI PER L'USO</i>	<i>I</i>	<i>49</i>
<i>DRIFTSANVISNING</i>	<i>S</i>	<i>53</i>
<i>DRIFTSINSTRUKS</i>	<i>DK</i>	<i>57</i>
<i>KÄYTTÖOHJE</i>	<i>N</i>	<i>61</i>
<i>DRIFTSVEJLEDNING</i>	<i>SF</i>	<i>65</i>
<i>ΟΔΗΓΙΕΣ ΛΕΙΤΟΥΡΓΙΑΣ</i>	<i>GR</i>	<i>69</i>
<i>INSTRUÇÕES DE SERVIÇO</i>	<i>P</i>	<i>73</i>

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Contens

	Page
Safety instructions	17
Proper use	18
Unpacking	18
Useful information	18
Drive unit	18
Assembly of boom and stand	19
Working with the unit	19
Dispersion tools: General notes	19
Dispersion tool S 25 N -	20
Dispersion tool S 25 NK -	20
Dispersion tool S 25 KR -	21
Dispersion tool S 25 KV -	21
Dispersion tool S 25 KV -- IL	22
Dispersion tool S 25 N - 8G	23
Dispersion tool S 25 N - 10G	24
Dispersion tool S 25 N - 10G - VS	24
Maintenance and cleaning	25
Accessories	25
Applicable standards and regulations	25
Guarantee	25
Technical data drive unit	26
Technical data for dispersion tools and properties	27
Technical data for dispersion tools and sealing types	28
Spare parts Drive unit T25 basic	80
Spare parts diagram T25 basic	81

Safety instructions

• Observe all pertinent safety instructions and directives, as well as work-safety and accident-prevention regulations for laboratory use. The unit must not be operated without a dispersion tool.

Never run dispersion tools dry, as the gasket and bearings will be destroyed if the tools are not cooled by the medium. The gasket and bearings are made of PTFE and rustproof steel; the following points should therefore be noted:

Chemical reactions of PTFE occur in contact with molten or solute alkali metals and alkaline earth metals, as well as with fine powders of metals in groups 2 and 3 of the periodic system at temperatures above 300 °C - 400 °C. Only elementary fluorine, chlorotrifluoride and alkali metals attack it; halogenated hydrocarbons have a reversible swelling effect.

(Source: Römpps Chemie-Lexikon and "Ulmann", Volume 19)

- Vibrations may loosen the hand screw on the flange of the drive unit. For your own safety check that the hand screw is secure and tighten it if necessary.
- Glass vessels must always be secured with a clamp to prevent them spinning. When working with flasks elastic couplers must be used to prevent fracture of the glass.
- Switch on the drive only at the lowest speed setting. Then gradually increase the speed to the desired level.
- Only dispersion and agitating tools approved by IKA may be used!
- Even in the case of repair work, the unit must only be opened by a qualified engineer. The mains plug must be pulled before opening the unit. The voltage-carrying components inside the unit may continue to carry a voltage for some time after the mains plug is pulled.
- Caution: Pull the mains plug before changing the dispersion heads!
- Caution: The integral switch **(2)** is only single-pole and is thus not suitable for safely isolating the unit from the supply!
- The ventilation slots on the drive must not be obstructed!
- The device does not start up again after a current failure.
- The T 25 basic may not be used to process inflammable or explosive media.

Proper use

The **ULTRA-TURRAX® T 25 basic** are dispersion units which, in connection with a dispersion tool, manufacture emulsions and dispersions. The dispersion devices should be operated on a laboratory stand.

Unpacking

The **ULTRA-TURRAX® T 25 basic** are entirely pre-assembled. Please unpack the device and/or dispersion tools carefully and inspect them for any damage. It is important that any damage is detected as soon as the items are unpacked. A record of the relevant facts may be required (railway, mail or carrier). The device delivery quantity includes: a driving unit, an extension arm, a hexagon socket screw, a tool set and an operating manual.

Useful information

You have acquired a high-quality product. Used in combination with an **S 25** series dispersion tool, the drive unit is a high-speed dispersing and emulsifying unit capable of handling free-flowing and liquid media in batches.

Dispersion is the dissolution and diffusion of a solid, liquid or gaseous phase in a continuum that is not consolute with that phase.


The rotor/stator principle:

Due to the high rotation speed of the rotor, the medium to be processed is automatically drawn axially into the dispersion

head and then forced radially through the slots in the rotor/stator arrangement. The high accelerations acting on the material produce extremely strong shear and thrust forces. In addition, high turbulence occurs in the shear gap between rotor and stator, which provides optimum mixing of the suspension.

The dispersion effectiveness is heavily dependent on the product of the shear gradient and the time the particles spend in the shear zone. The optimum range for the circumferential velocity of the rotor/stator arrangement is 10-24 m/s.

A processing time of a few minutes is usually sufficient to produce the desired fineness. Long processing times bring only insignificant improvements in the obtainable fineness; the energy expended serves merely to increase the temperature of the medium.

Drive unit

The open up a wide range of possibilities for dispersion technology in the conventional laboratory with a respective performance output of approximately 350 watts at 24,000 1/min in the case of the T25 basic.

The speed of the **T25 basic** dispersion unit set by using the differently adjusting wheel. In order to measure the speed of the dispersion unit an IKA **DZM control** tachometer may be connected (cf. Fig. 1) after removal of the cover.

Assembly of boom and stand

Overview (Figs. 1)

The extension arm is mounted as follows:

- insert extension arm **(1)** into flange;
- screw in the cylinder head screw **(2)**;
- firmly tighten cylinder head screw **(2)** using the SW4 right angled screwdriver.

Vibrations may loosen the cylinder head screw. Therefore, for your own safety, periodically check that the extension arm is secure. If necessary, retighten the cylinder head screw.

To ensure safe operation the driving units are fastened to the **R 1826 (8) stand base** by means of a **cross sleeve (7)**.

In order to enhance the stability of the mechanical arrangement, the driving units have to be mounted as close as possible to the stand pillar.

Working with the unit

Overview (Fig. 1)

Check that the voltage listed on the designation plate matches the mains voltage. Also pay attention to the ambient conditions listed in the technical data. The socket used must be fused to 16A.

The dispersion tool **(9)** is inserted into the driving unit until it stops and is then fixed by using the hand screw **(10)**. The

distance between the dispersion tool and the vessel bottom should not be less than 10 mm. The filling level of the medium may not be less than about 55 mm.

If these conditions are met the unit is ready for service when the mains plug has been plugged in. If this is not the case, safe operation cannot be guaranteed and the unit may be damaged.

The entire unit may also be arranged somewhat off centre in order to avoid any undesired air admission caused by the strong rotational turbulence.

The adjustment wheel **(4)** is set at the lowest speed before the unit is switched on. The speed which is set may be read by comparing the colors on the scale **(4)** and the adjustment wheel.

By using the adjustment wheel **(4)** the speed may be infinitely varied to meet the requirements of the medium to be processed. The no-load speed of the T25 driving unit may be set between 6,500 and 24,000 1/min.

The driving unit is switched on by means of the sliding switch **(6)**.

Dispersion tools: General notes

Overview (Fig. 3), (Fig. 4), (Fig. 5)

A dispersion tool consists of a shaft and a generator (dispersion head)

Shaft:

The shaft mainly consists of the shaft tube, the rotary shaft, the bearing and the shaft seal. There are three types of shaft for the **T25basic**:

S 25 N - / S 25 NK - / S 25 KR - / S 25 KV - / S 25 KV...-IL

The shaft cannot be used alone. Only combinations of shaft and generator can be used for working with the medium. The shaft is supplied with a cap (remove before use). Use this cap to protect shaft tubes from damage when not in service.

Generator:

The generator (dispersion head) consists of a rotor and a stator. These parts can be unscrewed from the shaft. Note that rotors and stators form matching pairs and must not be confused. Nine different dispersion heads are available for various applications:

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The dispersion tools you have ordered (e.g. S 25 NK - 19G) are always supplied fully assembled.

Never run dispersion tools dry, as the gasket and bearings will be destroyed if the tools are not cooled by the medium.

Dispersion tool S 25 N -

Overview (Fig. 5)

The bearing of the S 25 N - shaft is a friction bearing. This friction bearing is made of PTFE. The shaft has a dispersion head added when it leaves the factory.

The dispersion tool can be disassembled for cleaning as described below:

- Unscrew the rotor (**9**) from the shaft (**13**) using the generator and shaft key included in the tool set. The rotor is held fast by means of the generator key and screwed off from above by turning the shaft using the shaft key.
- Unscrew the stator (**10**) from the shaft (**4**) (**left-hand thread**). Insert flat key into the front-end slot of the shaft tube. By using the generator key the stator may now be screwed off in a clockwise direction.
- Pull the shaft (**13**) out of the shaft tube in a downward direction. The individual parts threaded on the shaft (**5, 6** and

7) may simply be pulled off.

- The slotted bearing (**7**) shall be replaced in the case of wear.

Reassemble the dispersion tool in the reverse sequence. When assembling ensure that the bearing bushing and the rotor are correctly installed.

Spare parts : Item	Designation
4	shaft tube
5	bushing
6	friction ring
7	slotted bearing
9	rotor
10	stator
13	shaft
25	protective cap

Sterilization :

After cleaning the dispersion tools the methods **A**), **B**) and **C**) are suitable for this dispersion tool.

(See maintenance, cleaning and sterilization of dispersion tools, page 25)

Dispersion tool S 25 NK -

Overview (Fig. 5)

The bearing of the **S 25 NK** shaft consists of a slide bearing and a ball bearing. The slide bearing is made of PTFE. The shaft is provided with a rotor when it leaves the factory.

The dispersion tool can be disassembled for cleaning as described below:

- Unscrew the rotor (**7**) from the shaft (**2**) using the generator and shaft key included in the tool set. The rotor is held fast by means of the generator key and screwed off from above by turning the shaft using the shaft key.

Risk of injury !

- Pull the shaft **(2)** out of the shaft tube in an upward direction. Push the bearing **(6)** out of the shaft tube **(1)** by using the bar included with the tool set.
- The bearing **(6)** shall be replaced in the case of wear.

Reassemble the dispersion tool in the reverse sequence.

spare parts :	Item	Designation
	1	shaft tube
	2	shaft
	6	bearing
	7	rotor
	15	protective cap

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After cleaning the dispersion tools the method **A** is suitable for this dispersion tool.

(See maintenance, cleaning and sterilization of dispersion tools, page 25)

Dispersion tool S 25 KR -

Overview (Fig. 5)

The bearing of the **S 25 KR** shaft consists of a ball bearing. The dispersion tool is sealed by means of rotary shaft seals **(5)** and an O-ring **(6)**. The shaft is provided with a dispersion head when it leaves the factory.

The dispersion tool can be disassembled for cleaning as described below:

- Unscrew the rotor **(10)** from the shaft **(14)** using the generator and shaft key included in the tool set. The rotor is held fast by means of the generator key and screwed off from above by turning the shaft using the shaft key.
- Unscrew the stator **(9)** from the shaft tube **(11)** (left-hand thread). Insert flat key into the front-end slot of the shaft tube. By using the generator key the stator may now be screwed off in a clockwise

direction.

- Pull the shaft **(14)** out of the shaft tube in an upward direction. The individual parts **(4, 6 and 12)** may now be simply removed from the shaft tube **(11)**.
- The rotary shaft seals **(5)** in the stator and the O-ring **(6)** must be replaced in the case of wear.

Reassemble the dispersion tool in the reverse sequence. In order to protect the rotary shaft seals push the bushing included in the tool set over the shaft thread. The guide angle of the bushing provides for easy assembly and prevents damage.

spare parts :	Item	Designation
	4	deep-groove ball bushing
	5	rotary shaft seal
	6	O-ring
	9	stator
	10	rotor
	11	shaft tube
	12	distance bushing
	14	rotary shaft assembly
	25	protective cap

Sterilization :

After cleaning the dispersion tools the method **A** is suitable for this dispersion tool. The methods **B**) and **C**) must not be used otherwise the ball bearing and seals will be destroyed.

(See Maintenance, cleaning and sterilization of dispersion tools, pages 25).

Dispersion tool S 25 KV -

Overview (Fig. 5)

The bearing of the **S 25 KV** shaft consists of two ball bearings and single-thrust slide ring seal. The slide rings are made of SiC. The shaft is provided with a dispersion head when it leaves the factory.

The dispersion tool can be disassembled for cleaning as described below:

- Unscrew the rotor **(9)** from the shaft **(2001)** using the generator and shaft key included in the tool set. The rotor is held fast by means of the generator key and screwed off from above by turning the shaft using the shaft key.
- Unscrew the stator **(8)** from the shaft **(2) (left-hand thread)**. Insert transverse rib of the shaft key into the front-end slot of the shaft tube. By using the generator key the stator may now be screwed off in a clockwise direction. In doing so, the dispersion tool and the stator should be held up in vertical position, thus preventing the spring of the slide ring sealing from falling out.
- In this position the parts of the slide ring seal can be removed from the shaft **(2001)**.
- Tap the shaft **(2001)** and the ball bearing **(4)** downward out of the shaft tube.
- The slide rings and O-rings must be replaced in the case of wear. The sliding surfaces of the slide and counter ring must be clean.

Reassemble the dispersion tool in the reverse sequence.

spare parts :	Item	Designation
	1	cylindrical pin
	4	deep-groove ball bushing
	6	Seeger circlip ring (crescent-shaped)
	7	porotective cap
	8	stator
	9	rotor
	10	O-ring 4,7 x 1,42
	11	O-ring 8,1 x 1,6
	12	O-ring 9,25 x 1,78
	14	deep-groove ball bushing
	15	ring
	25	protective cap
	2001	rotary shaft assembly
	2002	Friction ring seal assambly
	2003	shaft tube assembly

3005	slide ring
3006	counter ring
3003	pressure spring, right
3004	thrust collar
3008	counter fit mounting, compl

Make certain that the O-ring (10) and the thrust collar (3004) (Fig. 5) are properly fitted during assembly.

Sterilization :

After cleaning the dispersion tools the method **A**) is suitable for this dispersion tool. The methods **B**) and **C**) must not be used otherwise the ball bearing and seals will be destroyed.

(See Maintenance, cleaning and sterilization of dispersion tools, pages 25).

Dispersion tool S 25 KV - .. - IL

Overview (Fig. 5)

The bearing of the **S 25 KV - IL** shaft consists of a slide bearing, ball bearing and single-thrust slide ring seal. The slide rings are made of SiC. The shaft is provided with a dispersion head when it leaves the factory.

The dispersion tool can be disassembled for cleaning as described below:

- Unscrew the rotor **(9)** from the shaft **(2001)** using the generator and shaft key included in the tool set. The rotor is held fast by means of the generator key and screwed off from above by turning the shaft using the shaft key.
- Unscrew the stator **(8)** from the shaft **(2) (left-hand thread)**. Insert transverse rib of the shaft key into the front-end slot of the shaft tube. By using the generator key the stator may now be screwed off in a clockwise direction. In doing so, the dispersion tool and the stator should be held up in vertical position, thus preventing the spring of the slide ring sealing from falling out.

Risk of injury !

- In this position the parts of the slide ring seal can be removed from the shaft (**2001**).
- The slide rings and O-rings must be replaced in the case of wear. The sliding surfaces of the slide and counter ring must be clean.

Reassemble the dispersion tool in the reverse sequence.

spare parts :	Item	Designation
	1	cylindrical pin
	2	shaft tube
	4	deep-groove ball bushing
	6	Seeger circlip ring (crescent-shaped)
	7	protective cap
	8	stator
	9	rotor
	10	O-ring 4,7 x 1,42
	11	O-ring 8,1 x 1,6
	12	O-ring 9,25 x 1,78
	13	bearing
	25	protective cap
	2001	rotary shaft assembly
	2002	Friction ring seal assembly
	3005	slide ring
	3006	counter ring
	3003	pressure spring, right
	3004	thrust collar
	3008	counter fit mounting, compl

Make certain that the O-ring (10) and the thrust collar (3004) (**Fig. 5**) are properly fitted during assembly.

Sterilization :

After cleaning the dispersion tools the method **A**) is suitable for this dispersion tool. The methods **B**) and **C**) must not be used otherwise the ball bearing and seals will be destroyed.

(See Maintenance, cleaning and sterilization of dispersion tools, pages 25).

Dispersion tool S 25 N - 8G

Overview (Fig. 4)

The bearing of the **S 25 N - 8G** shaft is designed as a slide bearing. The slide bearing is made of PTFE.

The dispersion tool can be disassembled for cleaning as described below:

- Unscrew the rotor (**6**) from the shaft (**5**) using the flat and shaft key included in the tool set. The rotor is held fast by means of the flat key and screwed off from above by turning the shaft using the shaft key.
- Pull the shaft (**5**) out of the shaft tube in an upward direction.
- Take the grooved piston ring bearing (**7**) out of the shaft tube. The grooved piston ring bearing (**7**) and the check plate (**8**) must be replaced in the case of wear.

Reassemble the dispersion tool in the reverse sequence.

spare parts :	Item	Designation
	4	shaft tube
	5	rotary shaft
	6	rotor
	7	grooved piston ring bearing
	8	check plate
	25	protective cap

Sterilization :

After cleaning the dispersion tools the method **A**), **B**) and **C**) is suitable for this dispersion tool.

(See Maintenance, cleaning and sterilization of dispersion tools, pages 25).

Dispersion tool S 25 N - 10G

Overview (Fig. 4)

The bearing of the **S 25 N - 10G** shaft is designed as a slide bearing. The slide bearing is made of PTFE.

The dispersion tool can be disassembled for cleaning as described below:

- Unscrew the rotor **(6)** from the shaft **(10)** using the flat and shaft key included in the tool set. The rotor is held fast by means of the flat key and screwed off from above by turning the shaft using the shaft key.
- Pull the shaft **(10)** out of the shaft tube in an upward direction.
- Take the grooved piston ring bearing **(5)** out of the shaft tube. The grooved piston ring bearing **(5)** and the check plate **(4)** must be replaced in the case of wear.

Reassemble the dispersion tool in the reverse sequence.

spare parts :	Item	Designation
	4	check plate
	5	bearing bushing
	6	rotor
	7	centrifugal separator
	10	rotary shaft
	12	shaft tube
	21	protective cap
	25	protective cap

Sterilization :

After cleaning the dispersion tools the method **A)**, **B)** and **C)** is suitable for this dispersion tool.

(See Maintenance, cleaning and sterilization of dispersion tools, pages 25).

Dispersion tool S 25 N - 10G - VS

Overview (Fig. 4)

The bearing of the **S 25 N - 10G - VS** shaft consists of a ball bearing and a slide bearing. The slide bearing is made of PTFE.

The dispersion tool can be disassembled for cleaning as described below:

- Unscrew the rotor **(11)** from the shaft **(8)** using the flat and shaft key included in the tool set. The rotor is held fast by means of the flat key and screwed off from above by turning the shaft using the shaft key.
- Pull the shaft **(8)** out of the shaft tube in an upward direction.
- Take the bearing bushing **(10)** out of the shaft tube. The bearing bushing **(10)** has to be replaced in the case of wear.

Reassemble the dispersion tool in the reverse sequence.

spare parts :	Item	Designation
	5	O-ring
	6	deep-groove ball bushing
	8	rotary shaft
	9	Schleuderbuchse
	10	bearing bushing
	11	rotor
	13	shaft tube
	16	protective cap
	25	protective cap

Sterilization :

After cleaning the dispersion tools the method **A)** is suitable for this dispersion tool. The methods **B)** and **C)** must not be used otherwise the ball bearing will be destroyed.

(See Maintenance, cleaning and sterilization of dispersion tools, pages 25).

Maintenance and cleaning

Drive unit: The drive unit is maintenance-free but not immune from wear. The carbon brushes of the motor wear down over time and produce grooves in the collector. To clean the drive use only water with a detergent that contains tensides, or use isopropylalcohol for stubborn soiling.

Dispersion tools: The seals in the dispersion tools must be constantly monitored. In the event of leakage, the suction effect of the rotating shaft can cause the medium to penetrate as far as the drive unit. If liquid emerges from the side holes at the top of the shaft tube, stop work immediately and check the seals. The function of the dispersion heads depends on the condition of the sharp edges on the rotor and stator. These edges may be blunted very quickly in abrasive media, reducing the effectiveness of dispersion.

The dispersion tool is cleaned by operating it in a solvent which dissolves substance residues and is not harmful to the seal.

Because of the high flow speed, this is usually sufficient to clean the rotor and stator. The dispersion tool must be removed and cleaned immediately after use to prevent residues from adhering to the threads of the rotor and stator and to prevent the formation of undesirable bacteria cultures.

The dispersion tools may also be sterilized. The following methods may be used:

A) Chemical methods: Germicidal solutions (formalin, phenol, alcohol, etc.) are good general-purpose disinfectants. It is important that all disinfectant residues are then removed with sterilized water.. ***Make certain that the ball bearings, O-rings and shaft sealing rings are resistant to solvents.***

B) Sterilization with moist heat. This denotes the use of a steam jet pressurized to 2 bar at 120 °C.

C) Sterilization with hot air. Hot-air sterilization is usually performed at 160 °C or 190 °C. (30 min)

When ordering spare parts, please state the fabrication

number on the designation plate, the type of unit and the designation of the spare part.

Accessories

- **DZM control.o**
- **R 182** Cross sleeve
- **RH 3** Tensioner
- **R 1826** Plate stand
- **SV NS 29** Screw clamp
- **DK 25.11** Flow chamber

Applicable standards and regulations

EU Directives

EMC Directive: 89/336/EEC

Machines Directive: 89/392/EEC

EU standards

EN 61010-1:1993, VDE 0411-114

EN 50081:92, EN 55014:93, EN 50082:92, EN 60555:87

EN 292-1, EN 292-2, EN 60 204-1, EN 414

International standards

CAN/CSA C22.2 (1010.1:92), UL 3101-1

Guarantee

You have purchased an original IKA laboratory machine which meets the highest engineering and quality standards. In accordance with IKA guarantee conditions, the guarantee period is 24 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

Technical data

Speed range	rpm	6500	24000
Speed display		Scale	
Speed variation on load change	%	<6	
Perm. ambient temperature	°C	5 - 40	
Perm. humidity	%	80	
Perm. on time (drive unit)	%	100	
Protection to DIN 40050		IP20	
Power consumption	W	500	
Power output	W	300	
Design voltage	VAC	230 ±10%	
	or	VAC	115 ±10%
Frequency	Hz	50/60	
Noise level(wituout disp. tool)	dbA	73	
Dimensions (drive) WxDxH	mm	65 x 80 x 240	
Dimensions (boom)	mm	Ø13 / L160	
Weight	kg	1,6	
Contamination level		2	
Protection class		II	
Overvoltage category		II	
Operation at a terrestrial altitude	m	max. 2000 above sea level	

Technical data dispersing tools and properties

	Dispersion tool	Shaft tube-Ø mm	Ø Stator mm	Ø Rotor mm	Rotor / stator distance in mm	Peripheral speed m/sec	Recommended volume range ¹⁾ ml ²⁾	Available dispersion µm	Final emulsion fineness µm	Supplied sealing type S 25...	Comments	
	S 25 N – 8 G	8	8	6,1	0,25	7,7	0,3–5	10–50	1–10	N	with NS 14.5 metal polish	
	S 25 N – 10 G	10	10	7,5	0,35	9,4	1–50	10–50	1–10	N		
	S 25 N – 10 G-VS	10	10	7,5	0,35	9,4	1–50	10–50	1–10	N		
	S 25.. – 18 G	18	18	12,7	0,3	15,9	10–5000	10–50	1–10	N, KR, KV		
	S 25-NK – 19 G	19	19	12,7	0,3	16,0	25–2000	5–25	1–50	N		
	S 25... – 25 G	18	25	17	0,5	21,3	50–2500	15–50	1–10	N, KR, KV		
	S 25 KV. – 25 G	18	25	17	0,5	21,3	50–2500	15–50	1–10	KV		
	S 25.. – 25 F	18	25	18	0,5	22,6	100–2000	5–25	1–50	N, KR, KV		
	S 25 KV – 25 G-IL	18	25	17	0,5	21,3	Inline	10–30	1–10	KV		' UTL 25 basic
	S 25 KV – 25 F-IL	18	25	18	0,5	22,6	Inline	5-25	1–5	KV		' UTL 25 basic
	Additional rotors											
	– SR 18	–	–	12,7	0,3	15,9	10-1500	10 - 30	1 - 10		for unraveling	
	– SR 25	–	–	17	0,5	21,3	50–2500	10–30	1–10		for comminution	
	– RT 25	–	–	17	0,5	21,3	100–2000	5–25	1–50		for emulsions	
	– SW 18	–	–	12,8	0,35	15,9	10–500				f. viscous mater.	
	– SW 25	–	–	(G) 17 (F) 17	0,25 0,75	21,3	50–2500				f. viscous mater.	
	Assembly tools			Included in delivery quantity								

¹⁾ Insofar as flow is around seal

²⁾ Applies to 200 mPa · s

Technical data dispersing tools and sealing types

Sealing type	Maximum immersion depth in mm	Vacuum pressure in mbar	Material in contact with product	Suitable for use with solvent	pH range	Abrasive substances	Use up to °C ²	Capable of being sterilized	Recommended fields and industries					
									Bio ³	Ph ³	Med ³	Food ³	Ch ³	
S 25 N - 8 G	85	—	PTFE, 1.4571	yes	2–13	yes	+180°	yes all methods yes, wet chemic.	X	X	X	X	X	
S 25 N - 10 G	85	—												
S 25 N - 18 G	165	—												
S 25 N - 25 ..	165	—												
S 25 N - 10 G-VS	85	30 mbar												
S 25-NK - 19 G	165	—	PTFE, 1.4571	yes	2-13	yes	+120°	yes, wet chemic.	X	X	X	X	X	
S 25 KR	- 18 ..	185	FKM, 1.4571	no	2–13	no	+80°	yes, wet chemic.				X		
	- 25 ..	185												50 mbar
S 25 KV	- 18 ..	225	1.4571 FFPM/SIC 1.4571 " 1.4571 "	yes	2–13	(yes)	+220°	yes, wet chemic. yes, wet chemic. yes, wet chemic.	X	X		X	X	
	- 25 ..	225												1 mbar
	- 25 ... -IL	85												to 6 bar

¹Abrasion may be expected as of hardness 6; ²applies to medium temperature, without adding warming generated by shear.

³Bio = biology, biotechnology; Ph = pharmaceutics, cosmetics; Med = medical diagnostics; Food = foodstuffs; Ch = chemistry, process engineering.

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