Bedienungsanleitung Operating Instructions Mode d'emploi

SONOREX SUPER und SONOREX LONGLIFE

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SONOREX SUPER and SONOREX LONGLIFE

ultrasonic cleaners for aqueous cleaning agents

Cleaning units, accessories and cleaning agents have to be strictly used according to the instruction manual or the corresponding product information.

Keep cleaning units and cleaning agents away from children and from persons not instructed for correct use by means of these instructions.

In case of improper use, the manufacturer does not take any responsibility for safety and proper functioning.

Unauthorized changes/modifications will result in expiration of the CE-conformity.

If service is required, please contact your dealer or the address mentioned below.

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Contents

1	Strictly observe the following general instructions for safety and correct operation!	16
2	Purpose	17
2.1	Preparation of medical instruments	17
2.2	EC-conformity	17
3	Start	17
3.1	Degassing	17
3.2	Cleaning	18
3.3	Ultrasound ON / OFF	19
3.4	Units with heating	19
3.5	Emptying the ultrasonic tank	20
4	Recommendations	20
4.1	Cleaning agents	20
4.2	Cleaning and care	21
5	Disfunctions.....	22
6	Dosing table for SONOREX – compact ultrasonic cleaners	23
7	Technical data	24 and 39
8	Electromagnetic environmental conditions.....	25

Copyright & limit of liability

This document may not be reproduced, either in full or in extracts, without prior approval of BANDELIN electronic GmbH & Co. KG, hereinafter referred to as BANDELIN.

The German original is the binding version of this document. Any difference in the translation is not binding and has no legal effect. In case of any discrepancies between the translation and the original, the original version has priority.

BANDELIN accepts no responsibility and liability for damage caused by improper handling or usage contrary to the intended purpose.

The documentation was created with great care. Liability for indirect and direct damage arising because of incomplete or erroneous information in this documentation as well as in delivery and usage is excluded.

© 2009

BANDELIN *electronic* GmbH & Co. KG

Strictly observe the following general instructions for safety and correct operation!

- Keep cleaning units and cleaning agents away from children and also away from persons not instructed for correct use by means of these instructions.
- Do not immerse parts of the body (hand, feet) nor creatures or plants into the ultrasonic tank. During ultrasonic cleaning do not dip your hands into the cleaning liquid.
- Avoid any contact of your eyes and/or your skin with the cleaning liquid.
- Do not ingest nor breathe in the cleaning liquids.
- Never run the unit without liquid!
- Do not fill hot water into the ultrasonic tank. The filling temperature must not exceed 50 °C.
- Do not put parts to be cleaned into the tank with aggressive soiling such as acids, chloride ions or similar substances!
- Do not fill the stainless steel oscillating tanks with aggressive cleaning agents such as acids or chemicals which contain or release chloride ions (some disinfectants, dish-washing detergents, household cleaners, saline solutions).
- Do not fill the stainless steel oscillating tanks with combustible liquids.
- Use aggressive cleaning agents only in inset beakers or non-perforated plastic insert tanks.
- After removing the packing, check the unit's casing and control elements for any signs of damage. In such case, do not connect the unit to the electrical power supply, but report the damage immediately and in writing to the carrier and your supplier. Save the original packing.
- Place the unit in a horizontal position on a dry, solid base.
- Plug the unit only to a socket with grounded safety contact.
- Keep the unit's surface and control elements clean and dry!
- Replace defective parts only with original **SONOREX** parts!
- During preheating, it is necessary to stir the cleaning liquid at least every 15 minutes or switch on ultrasound (retardation of boiling – danger of scalding)!
- Ultrasonic cleaning produces cavitation noise. If working within distances of 2 meters from the unit for longer periods of time, be sure to wear ear plugs to avoid damage to your hearing!
- Do not place anything on the bottom of the oscillating tank. Use appropriate **SONOREX** accessories (insert basket, positioning lid with insert vessel, etc.).
- The ultrasonic energy in the oscillating tank heats up the cleaning liquid without additional heating! Check the temperature regularly when cleaning temperature-sensitive parts.
- Do not move or carry units filled with liquid.
- Do not allow the unit to run without supervision.

Principle of ultrasonic cleaning

PZT oscillators underneath the tank convert electrical energy into mechanical vibrations.

SONOREX units produce oscillations of 35 kHz in the cleaning liquid. Tiniest vacuum bubbles are formed which then implode (cavitation). This cavitation principle removes dirt and grease from the deepest pores - "electronic scrubbing" - at all points that come into contact with the cleaning liquid.

2 Purpose

2.1 Preparation of medical instruments

- In this case, the ultrasonic bath is used as an ultrasonic cleaning device in conjunction with disinfection and/or cleaning solutions, in order to speed up the disinfection and/or cleaning process.
- The instruments must not be placed on the tank bottom. The basket prevents damage to the instruments and the floor of the bath. Ensure that the instruments are spread out in the bath. Overloading the basket will reduce the cleaning efficacy. Pliers or scissors should be completely opened, or disassembled if possible. Ensure that the instruments are completely submerged in the solution. Air from cavities must be able to escape.

2.2 EC-conformity

Our ultrasonic cleaning units are declared as „Medical products“, and are in compliance with the EC-designation criteria of the European guidelines such as:

- guideline „Medical products“
- „EMC-guideline“
- guideline „Electromagnetic compatibility“

in their current versions.

The EC-conformity declaration can be requested from the manufacturer.

* Exclusion: refer to the general units survey on page 39.

3 Start

- Make sure that the ball valve of the unit is closed!
- Fill the oscillating tank up to the filling mark.
- Connect the unit to a grounded socket.

3.1 Degassing

Degassing the cleaning liquid enhances the cleaning effect.


Cleaning and contact liquids contain dissolved gases (for example, oxygen). Sonicate cleaning and contact liquid, which have just been poured into the oscillating tank or have been left in it for a longer period of time, for approx. 5 to 15 minutes before starting the cleaning process. Degas oscillating tanks with a volume of more than 10 l for roughly 30 minutes. During the degassing process the cavitation noise changes. As the loud degassing noise decreases towards the end of the process, the unit appears to run more quietly. Lower noise levels indicate the end of the degassing process and an improved cleaning impact. They do not mean a reduction of ultrasonic energy.

3.2 Cleaning

- Be sure that the objects in the tank are completely covered with liquid.
- Normally the objects may be cleaned **directly** in the ultrasonic tank with an insert basket.
- In the following cases clean objects **indirectly** in inset beakers or non-perforated insert tanks to protect the stainless steel oscillating tank:
 - when chemically aggressive cleaning liquids (for example acids, etc.) are used.
 - when chemically aggressive substances (for example when cleaning racks from film developing machines) or
 - when abrasive substances (for example polishing paste, quartz, sand) are to be removed.

DIRECTLY in the oscillating tank:


- Fill the oscillating tank to 2/3 with cleaning liquid (TICKOPUR, STAMMOPUR see dosage on label). Improperly filled units can suffer damage (see point 5 – dosing table).
- Use always an insert basket or utensil holder.
- Parts must not be placed on the tank bottom (Exemption basket K6 and SH 7).

INDIRECTLY in inset beakers or non-perforated insert plastic tub


⇒ **SONOREX** Accessories

- Fill the oscillating tank with a tenside-containing contact liquid (for example TICKOPUR R 33 or STAMMOPUR RD 5 diluted with water) to achieve a good transmission of ultrasound.
- Fill the inset beakers or insert plastic tub with cleaning liquid (for example TICKOPUR, STAMMOPUR).
- When sonicating small amounts of combustible liquids in insert beakers, the country-specific guidelines/regulations in the particular valid versions have to be observed.
- Suspend inset beakers in the positioning cover, do not place them on the tank bottom.
- Inset beakers or insert plastic tub must be immersed in the contact liquid to a depth of at least 2 cm. See to the tank filling to 2/3.
- Several inset beakers with different cleaning liquids can be sonicated at the same time.

When inserting the objects to be cleaned, be sure to observe the following:

- Place the dirtier side towards the bottom.
- Do not stack - ultrasonic energy will be absorbed.
- Sensitive pieces may not be allowed to touch one another. Avoid damage by using plastic insets, plastic insert baskets or the knobbed silicone mats in insert basket.
- Be sure that there are no air bubbles in hollow spaces (for example pocket holes).
- Check the temperature when cleaning heat-sensitive objects.

3.3 Ultrasound ON / OFF

Cleaning time


Cleaning time must be kept as short as possible. Depending on the degree of soiling and the cleaning liquid used (see the data on its label), the cleaning time varies between 1 and 15 minutes.

Continuous operation (∞ position)

- Turn the knob clockwise until it clicks into ∞ position → the green pilot lamp lights up.
- The unit does not switch off automatically, turn the knob back to "0" in order to switch off.

Short time operation

- Turn the knob clockwise, beyond continuous operation, to the required cleaning time → the green pilot lamp lights up.
- The timer switches off automatically after run-off of time.
- You can reduce the time or switch the unit off by turning the knob back.

Leads for timer:

- In the switched off state, the unit may remain plugged in. To disconnect it, pull the mains plug.
- The clock must not be turned to the left in the switched off state ("0" position)!

Important leads for temperature:

- Warmed-up liquids improve the cleaning efficiency.
- However, the effect of the ultrasonic cavitation decreases at higher temperatures¹. Experience has shown that best cleaning results are achieved by using an aqueous solution with a bath temperature of 50 to 60 °C.
- Ultrasonic energy also heats the cleaning agent. Continuous irradiation and/or the covering of the oscillating tank may cause the temperature of the liquid to rise more quickly and exceed the thermostat's default setting (80 °C).
- **Do not** heat disinfectants.
- When reaching a critical temperature, the ultrasonic power is reduced to protect the electronic components inside the ultrasonic device. This way, the internal temperature is prevented from further rising.


3.4 Units with heating

Units marked with „H” have a built-in heating.


The heating system runs independently from the ultrasonic energy.

- Tip switch operation (fixed temperature of 65 °C): Heating switch on – control lamp is on. When the bath temperature of approx. 65 °C is reached, the heating will switch off, the control lamp remains on.
- Set the temperature knob clockwise to the selected temperature → the thermostatic regulation is in operation:
 - the yellow and white pilot lamps light up.
 - the yellow pilot lamp extinguishes when the selected temperature is reached.

RK 31 H: Tip switch

RK 52 H/CH and upwards: temperature knob

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
Units with drain


Ball valve

Turn the lever to the outflow direction.

Units without drain

Pour out the liquid over the rear left tank corner. Avoid splashing from below and into the socket.

4 Recommendations

4.1 Cleaning agents

- Only use approved and especially developed cleaning agents, e.g. TICKOPUR, TICKOMED, STAMMOPUR in order to obtain a high cleaning efficiency. Tap water without any additives does not clean sufficiently.
- When using the cleaning agents you must always strictly observe the safety instructions in the product information.
- Keep cleaning agents away from children and from persons not instructed for correct use by means of the product information.
- Do neither ingest nor inhale the cleaning agents. Avoid any contact of your eyes and/or skin with the cleaning agents.
- Replace used liquids and disinfectants, do not “freshen them up” with an added dose.
- Layers of dirt on the tank bottom or in the insets (abrasives and polishing agents, for example, as well as other deposits) reduce the cleaning efficiency. Empty the oscillating tank and clean it regularly. Replace the cleaning liquid.
- The cleaning and contact liquids in the ultrasonic oscillating tank may not exceed a maximum operating temperature of 100 °C.
- Distilled or deionized water without additives is unsuitable for ultrasonic cleaning directly in the oscillating tank. They lack the necessary cavitation properties to form a uniform ultrasonic field. Irregular cavitation can damage the oscillating tank or the objects being cleaned. (You may use distilled or deionized water in the inset beakers and non-perforated insert tubs.)
- Be careful when working with aggressive cleaning agents in insets or non-perforated insert tanks; avoid splashing into the contact liquid or onto stainless steel surfaces; if necessary, replace the contact liquid immediately, clean the surfaces and wipe them dry.

The manufacturer does not take any responsibility for damage of equipment or damages at the parts to be cleaned caused by inappropriate cleaning chemicals.

4.2 Cleaning and care

General

- Be sure to pull the plug before cleaning the unit.
- Do not shower the units or immerse them in water; instead, wipe them with a damp cloth, in- and outside.
- Rinse out the stainless steel tank frequently and wipe it dry.
- Do not use steel wool, scrapers or graters.
- If rings should appear in the stainless steel tank after a longer period of operation, remove them with a customary stainless steel cleaner without scouring additives.
- Metal parts remaining on the stainless steel surface, as well as particles of rust, can penetrate the passive protective layer of the stainless steel, thereby "activating" it and causing it to rust. This rust provokes pittings on the stainless steel surface. Therefore, remove all remaining metal parts such as screws, metal chips and the like; remove minor rust spots immediately with a soft cloth and customary stainless steel cleaners without scouring additives.

Medical use

- If contaminated medical instruments are being prepared in the ultrasonic bath, hygienic safety following their use is important. In order to avoid cross-contamination by the settlement of micro-organisms especially in the edge of the bath and the drain, but also on the work surfaces these areas should be cleaned and disinfected regularly with a VAH-certified surface disinfectant or another disinfectant found to be effective, according to the hygienic plan.

SONOREX units do not require maintenance.

5 Disfunctions

- Oscillations are too weak and non-uniform, the oscillation noise is too loud; the cleaning effect is insufficient:
 - Has the liquid been properly degassed? ⇒ sonicate for 15 minutes.
 - Is the unit overloaded? ⇒ remove some parts.
 - Non-uniform noise (wobbling) is not a defect ⇒ change the fluid level.
- Heating defect:
No objection to operate the unit without the heating.
- Minor erosion marks at the tank bottom? ⇒ common wear and tear effect.
The unit is faultless.
- The mains fuse can be released. ⇒ The unit is equipped with an overvoltage arrester which protects it from voltage peaks; the mains fuse can be released during overcharge.
- Repairs should only be carried out by authorised service.
It is imperatively essential to disconnect the mains plug before opening the device!
Danger of electric shock caused by live parts in the appliance!
The manufacturer will disclaim all liability in case of unauthorized intervention!
- Replace defective parts only by original **SONOREX** parts.

Please inform us about disfunctions in writing.

6 Dosing table for SONOREX – compact ultrasonic cleaners

You can order this dosing table (DIN A 4-size) free of charge or download it as PDF-file from the following internet address::

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Values have been rounded up mathematically.

number printed in normal type: water

number printed in **bold** type: concentrate

Fill the oscillating tank to 1/3 with water and add the appropriate quantity of concentrate. Then add water up to the filling mark and sonicate for approx. 15 minutes for degassing.

or

According to the dosing table, fill the oscillating tank with water, add the appropriate quantity of concentrate, mix the solution for a short while and sonicate for approx. 15 minutes for degassing.

Calculation of other volumes

Example:

- 10 litres solution ready for use
- 2,5 % dosage of the concentrate

$$\frac{10 \text{ l} \times 2,5 \%}{100 \%} = 0,25 \text{ l concentrate}$$

$$10 \text{ l} - 0,25 \text{ l} = 9,75 \text{ l water}$$

7 Technical data

SONOREX compact units are RFI-proofed and ζ ϵ marked.
Compliance with the limits according to EN 60601-...

Nominal voltage: 230V~ 50/60 Hz, (115 V on request)
cable length 2 m
High frequency: 35 kHz
Degree of protection: IP 32
Oscillating tank: stainless steel
Serial number: see type label on the rear side of the unit

Environmental conditions in accordance with EN 61 010-1 (IEC 1010-1)

Degree of soiling: 2 according to IEC 60664-1:11.02
Excess voltage category: II
Permissible ambient temperature: 5 - 40 °C
Permissible relative humidity up to 31 °C: 80 %
Permissible relative humidity up to 40 °C: 50 %
Dew not permitted.

Information for the use as a medical product

Denotation: ultrasonic cleaning device
UMDNS-nomenclature (ECRI): 14-263
Purpose determination: ultrasonic-disinfection and cleaning of medical and dental instruments
Classification (Medical Devices Directive 93/42/EEC, appendix IX): class I; active, not invasive, not implantable medical device
Type, model, serial number, year of manufacture: see type designation plate on rear of device

Information according to IEC 601-1:

Protection class: class I

8 Electromagnetic environmental conditions

Electromagnetic Emissions		
The SONOREX is intended for operation in the following electromagnetic environments. The user should make sure that the SONOREX is used in such an environment.		
Interference emission measurements	Conformity	Electromagnetic Environment – Guidelines
HF emissions as per CISPR 11	Group 1	The SONOREX uses HF energy exclusively for internal functioning. The HF emission is therefore very small and it is unlikely that adjacent electronic units are disrupted.
HF emissions as per CISPR 11	Class B	The SONOREX is suitable for use in facilities outside the residential area and similar places that are directly connected to a public power supply network which also supplies buildings used for residential purposes.
Emission of harmonics as per IEC 61000-3-2	Class A	
Emission of voltage fluctuations/ flickers as per IEC 61000-3-3	Coincides	

Electromagnetic Stability			
The SONOREX is intended for operation in the following electromagnetic environments. The user should make sure that the SONOREX is used in such an environment.			
Interference immunity tests	IEC 60601 test level	Consistency level	Electromagnetic Environment – Guidelines
Discharge of static electricity (DSE) as per IEC 61000-4-2	± 6kV contact discharge ± 8kV air discharge	± 6kV contact discharge ± 8kV air discharge	Floors should be made of wood or concrete or covered with ceramic tiles. If the floor is made of synthetic material, the relative humidity must at least be 30%.
Rapid transient electrical disturbance variables/bursts as per IEC 61000-4-4	± 2kV for power lines ± 1kV for inbound and outbound lines	± 2kV for power lines ± 1kV for inbound and outbound lines	The quality of the supply voltage should correspond to the typical (business or) hospital environment.
Surges as per IEC 61000-4-5	± 1kV differential mode voltage ± 2kV common-mode voltage	± 1kV differential mode voltage ± 2kV common-mode voltage	The quality of the supply voltage should correspond to the typical (business or) hospital environment.
Voltage dip, short-term interruptions and fluctuations of supply voltage as per IEC 61000-4-11	<5% U_T (>95% decline of U_T for 1/2 period) 40% U_T (60% decline of U_T for 5 period) 70% U_T (30% decline of U_T for 25 period) <5% U_T (>95% decline of U_T for 5s)	<5% U_T (>95% decline of U_T for 1/2 period) 40% U_T (60% decline of U_T for 5 period) 70% U_T (30% decline of U_T for 25 period) <5% U_T (>95% decline of U_T for 5s)	The quality of the supply voltage should correspond to the typical (business or) hospital environment. Continuation of functioning of the SONOREX is assured even after interruptions of the power supply occur.
Magnetic field with supply frequency (50/60Hz) as per IEC 61000-4-8	3A/m	3A/m	Magnetic fields at network frequency should correspond to the typical values found in (office and) hospital environments.
Observation: U_T is the network AC voltage before application of test level.			

Electromagnetic Stability

The SONOREX is intended for operation in the following electromagnetic environments. The user should make sure that the SONOREX is used in such an environment.

Interference immunity tests	IEC 60601 test level	Consistency level	Electromagnetic Environment – Guidelines
			Portable and mobile radios should not be used at a distance to the SONOREX including the lines closer than the recommended safety distance calculated according to the equation applicable to the transmitter frequency. Recommended safety distance:
Directed HF disturbance variables as per IEC 61000-4-6	$3 V_{\text{eff}}$ 150 KHz to 80 MHz	$V_1 = 3 V_{\text{eff}}$	$d = 1,2\sqrt{P}$
Irradiated HF disturbance variables as per IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	$E_1 = 3 \text{ V/m}$	$d = 1,2\sqrt{P}$ for 80 MHz to 800 MHz $d = 2,3\sqrt{P}$ for 800 MHz to 2,5 GHz
			With P as rated power of the transmitter in Watt (W) according to specifications of the transmitter manufacturer and d as the recommended safety distance in meters (m). The field strength of stationary radio transmitters should be less than the consistency level at all frequencies according to an on-site examination. Disturbances are possible in the environment of units which carry the following graphical symbol: 

Recommended safety distances to portable and mobile HF communication units

The SONOREX is intended for operation in an electromagnetic environment in which the irradiated HF disturbance variables are controlled. The user of the SONOREX can thus help prevent electromagnetic disturbances by maintaining the minimum distance between portable and mobile HF communication units (transmitters) and the SONOREX – depending on the output power of the communication unit:

Rated power of the transmitter in W	Safety distance according to transmitter frequency in m		
	150 kHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,37	0,37	0,74
1	1,17	1,17	2,33
10	3,70	3,70	7,37
100	11,70	11,70	23,33

For transmitters whose maximum rated power is not indicated in the table above, the recommended safety distance d in metres (m) can be calculated using the equation belonging to the respective column, where P is the maximum rated power of the transmitter in Watt (W) according to the specification of the transmitter manufacturer.

Observation 1: At 80 MHz and 800 MHz the higher value applies.

Observation 2: To calculate the safety distance of transmitters in the frequency range of 80 MHz to 2.5 GHz, an additional factor of 10/3 is used to reduce the probability that a mobile/portable communication unit will result in a disturbance.

Observation 3: These guidelines may not be applicable in all cases. The propagation of electromagnetic quantities is affected by absorption and reflections from buildings, objects and people.

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