## KERN **EMS 3000-2**

Entry level laboratory balance with large weighing plate





Measuring system	
Adjustment options:	External calibration
Eccentric loading at 1/3 [Max]:	50 mg
Linearity:	50 mg
Readability [d]:	10 mg
Recommended adjusting weight:	3 kg (F2)
Repeatability:	20 mg
Resolution:	300000
Stabilisation time under laboratory conditions:	3 s
Tare range:	3 kg
Warm up time:	2 h
Weighing capacity [Max]:	3 kg
Weighing system:	Strain gauge
Weighing units:	g oz dwt ozt

Display digit height (large):	2,500 cm
Counting	
Counting resolution:	150000
Mininum piece weight at piece counting (Laboratory:	20 mg
Mininum piece weight at piece counting (Normal):	200 mg

Display

Functions	
Counting function:	yes
Percentage determination:	yes

Environmental conditions	
Maximum humidity:	80 %
Maximum operating temperature:	35 °C
Minimum ambient temperature:	5°C
Minimum humidity:	40 %

1 Owel Supply	
Battery:	9 V Block
Input voltage:	220 V - 240 V AC 50 Hz
Mains adapter:	Mains adapter external
Mains adapter/adapter included:	EURO CH US UK
Operating time:	40 h
Rechargeable battery:	Rchrg. battery optional

Services (optional)	
DAkkS Certificate:	963-127

Category	
Brand:	KERN
Category:	Balances
Product Group:	School balance

Delivery:	24 h
Dimensions packaging (WxDxH):	310 x 234 x 164 mm
Gross weight:	1,980 kg
Net weight:	1,575 kg
Shipping method:	Parcel service

Plastic
200 x 280 x 48 mm
Plastic
175 x 190 mm

## Piktogramme

Packaging & shipping

## STANDARD



## **OPTION**

