

DATASHEET

For in vitro research only.
Not for diagnostic or clinical use.

Cat.No. 317400

Transparent Sealing Tape Storage 100 sheets

BIOVENDIS

Biovendis Products GmbH
Joseph-Meyer-Str . 13-15
68167 Mannheim

Tel. + 49 - 621 - 43 00 52 52
Fax + 49 - 621 - 43 00 52 50
E-Mail info@biovendis.de

PRODUCT

Transparent Sealing Tape Storage

ORDERING INFORMATION

Cat. No. 317400 - Transparent Sealing Tape Storage

PLATE MATERIAL

Polypropylene (PP)

SEALING CONDITIONS

Adhesive Sealing Tape.

SEAL INTEGRITY RANGE

-20°C to 80°C;

This product is not recommended for high temperature applications like PCR or storage at temperatures below -20°C!

DESCRIPTION

The clear adhesive sealing tape offers a convenient and cost effective adhesive cover for long and short-term storage.

Tape creates a good reseal on top of pierced tapes. End tab for easy seal application and removal

STORAGE

Seals should be stored in a cool, dry place. Avoid direct exposure to sunlight. It is recommended to use seals within one year from date of purchase.

SHEET SIZE

133.2 x 83.0 mm (complete)

120.0 x 83.0 mm (adhesive section)