

# **THERMOSPRINT®** **Transparent Sealing Tape Storage**



100 sheets

## **Product**

THERMOSPRINT Transparent Sealing Tape Storage

## **Ordering Information**

Cat. No. 317400 - Transparent Sealing Tape Storage

## **Plate Material**

Polypropylene (PP)

## **Sealing Conditions**

Adhesive Sealing Tape.

## **Seal Integrity Range**

-20°C to 80°C;

**This product is not recommended for high temperature applications like PCR or storage at temperatures below -20°C!**

## **Description**

The clear adhesive sealing tape offers a convenient and cost effective adhesive cover for long and short-term storage.

Tape creates a good reseal on top of pierced tapes.

End tab for easy seal application and removal

## **Storage**

Seals should be stored in a cool, dry place. Avoid direct exposure to sunlight. It is recommended to use seals within one year from date of purchase.

## **Sheet Size:**

133.2 x 83.0 mm (complete)

120.0 x 83.0 mm (adhesive section)