

**THERMOSPRINT®**



**Transparent Heat Sealing Tape extra strong**

100 sheets

**Product**

THERMOSPRINT Transparent Heat Sealing Tape extra strong

**Ordering Information**

Cat. No. 317700 - Transparent Heat Sealing Tape extra strong

**Plate Material**

Polypropylene (PP)

**Sealing Conditions**

2-3 seconds

1.5 – 2.0 seconds at 165°C to 170°C.

**Seal Integrity Range**

-80°C to 120°C;

(120°C with adequate pressure from cyclor lid)

**Description**

The tape is a clear polyester polypropylene laminate appropriate for fluorescence and colorimetric assay detection.

Good optical clarity suitable for PCR and QPCR applications.

Weald seal cannot be peeled or pierced. Due to its strong and unpeelable seal properties, it is also suitable for storage and disposal of hazardous material.

**DMSO Resistance**

Good solvent resistance. Recommended for use with DMSO (dimethyl sulphoxide)

**Storage**

Seals should be stored in a cool, dry place. Avoid direct exposure to sunlight. It is recommended to use seals within one year from date of purchase.