

What version of kit to choose

The BUGS'n BEADS™ kits are delivered in three different versions developed to target different groups of bacteria. The S-version binds with superior performance to bacteria commonly associated with food poisoning. The C-version binds with superior performance to species commonly associated with serious infections such as meningitis and bacteraemia. The U-version binds to bacteria such as Chlamydia, Mycobacterium and MRSA.

Bacteria binding to BUGS'n BEADS S-version (GK1) are:

Aeromonas
Bacillus*
Campylobacter
Citrobacter
Clostridium
E.coli
Enterobacter
Enterococcus
Hafnia
Klebsiella
Lactobacillus
Listeria
Morganella
Providencia
Proteus
Salmonella
Serratia*
Shewanella
Shigella
Vibrio
Yersinia

Bacteria binding to BUGS'n BEADS C-version (GK2) are:

Acinetobacter
Alcaligenes
Bacillus*
Bordetella
Borrelia
Chlamydia
Clostridium
Corynebacterium
E. coli*
Enterobacter
Haemophilus
Helicobacter
Klebsiella
Listeria
Micrococcus
Mycobacterium
Neisseria
Propionebacterium
Proteus
Pseudomonas
Salmonella*
Serratia
Streptococcus*
Staphylococcus*
Yersinia

Bacteria binding to BUGS'n BEADS U-version (GK3) are:

Bacillus
Chlamydia
Listeria
Mycobacterium
Mycoplasma
Neisseria
Staphylococcus
Streptococcus
Ureaplasma

* Bacteria will bind, but using another kit-version results in binding with higher affinity.